


Anti-Cnpase Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P09543

Item Number: CM0003321



Introduction

High-specificity rabbit monoclonal antibody targeting human CNPase (P09543), validated for WB, IF-P, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for myelin biology and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	CNPase Rabbit mAb
Remarks/Alias	CNP1; HLD20; CNPase
Species	Human
GeneID (Human)	1267
GeneID	1267
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 322-421 of human CNPase (P09543)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09543
Protein Weight	48kDa
Shipping	Ice bag