


Anti-Cnpase Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P09543

Item Number: CM0014807



Introduction

CamelBio offers a high-quality rabbit polyclonal antibody targeting human CNPase (also known as CNP1/CNP). Validated for WB, IF-P, and ELISA with cross-reactivity to mouse and rat, it supports myelin, neurodegeneration, and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	CNPase Rabbit pAb
Remarks/Alias	CNP1; HLD20; CNPase
Species	Human
GeneID (Human)	1267
GeneID	1267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-421 of human CNPase (NP_149124.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09543
Protein Weight	48kDa
Shipping	Ice bag