

Anti-Cnpy2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y2B0

Item Number: CM0012995

Introduction

Rabbit polyclonal antibody against human CNPY2 (MSAP, TMEM4, ZSIG9), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for neurite outgrowth research.

[Learn More](#)

Parameter	Value
Product Name	CNPY2 Rabbit pAb
Remarks/Alias	MSAP; TMEM4; ZSIG9; HP10390; CNPY2
Species	Human
Gene ID (Human)	10330
Gene ID	10330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of human CNPY2 (NP_055070.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2B0
Protein Weight	21kDa
Shipping	Ice bag