

Anti-Cnpy3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9BT09

Item Number: CM0025063

Introduction

Rabbit polyclonal antibody against human CNPY3 (PRAT4A), validated for WB, IHC-P, IP, and ELISA. Cross-reacts with mouse. Essential for innate immunity and ER protein folding studies.

[Learn More](#)

Parameter	Value
Product Name	CNPY3 Rabbit pAb
Remarks/Alias	CAG4A; DEE60; ERDA5; TNRC5; EIEE60; PRAT4A; CNPY3
Species	Human
GeneID (Human)	10695
GeneID	10695
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-278 of human CNPY3 (NP_006577.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BT09
Protein Weight	31kDa
Shipping	Ice bag