


Anti-Cnr1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P21554

Item Number: CM0016536



Introduction

Rabbit polyclonal antibody against human CNR1 (CB1), suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Useful for cannabinoid receptor research in neuroscience and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	CNR1 Rabbit pAb
Remarks/Alias	CB1; CNR; CB-R; CB1A; CB1R; CANN6; CB1K5; CNR1
Species	Human
GeneID (Human)	1268
GeneID	1268
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 401-472 of human CNR1 (NP_001153698.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21554
Protein Weight	53kDa
Shipping	Ice bag