

Anti-Cntfr Rabbit Polyclonal Antibody For Wb, Elisa - P26992

Item Number: CM0029828



Introduction

Rabbit polyclonal antibody targeting human CNTFR, suitable for WB and ELISA. Reacts with human, mouse, and rat. Enables study of CNTF signaling, motor neuron survival, and neuroprotection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	CNTFR Rabbit pAb
Remarks/Alias	CNTFR
Species	Human
GeneID	1271
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-342 of human CNTFR (NP_671693.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26992
Protein Weight	41kDa
Shipping	Ice bag