


Anti-Cntn1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12860

Item Number: CM0023601



Introduction

Rabbit polyclonal CNTN1 antibody, also known as F3/GP135, validated for WB, IHC-P, ELISA in human, mouse, and rat. Ideal for neural development research studies.

[Learn More](#)

Parameter	Value
Product Name	CNTN1 Rabbit pAb
Remarks/Alias	F3; GP135; MYPCN; CMYP12; CNTN1
Species	Human
GeneID (Human)	1272
GeneID	1272
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CNTN1 (NP_778203.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12860
Protein Weight	113kDa
Shipping	Ice bag