

Anti-Cntn3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9P232

Item Number: CM0016179



Introduction

Rabbit polyclonal antibody against human CNTN3 (Contactin-3, BIG-1, PANG) validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for neural cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	CNTN3 Rabbit pAb
Remarks/Alias	PCS; PANG; BIG-1; CNTN3
Species	Human
GeneID (Human)	5067
GeneID	5067
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human CNTN3 (NP_065923.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9P232
Protein Weight	113kDa
Shipping	Ice bag