

Anti-Cntn4 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Iwv2

Item Number: CM0014056

Introduction

CNTN4 rabbit polyclonal antibody validated for WB and ELISA. Targets human Contactin-4 (BIG-2), a neural cell adhesion molecule involved in neurite outgrowth and synaptogenesis. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CNTN4 Rabbit pAb
Remarks/Alias	AXCAM; BIG-2; CNTN4
Species	Human
GeneID (Human)	152330
GeneID	152330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 448-697 of human CNTN4 (NP_001193885.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8Iwv2
Protein Weight	113kDa
Shipping	Ice bag