


Anti-Coagulation Factor Iii/Tissue Factor Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13726

Item Number: CM0018174



Introduction

Rabbit polyclonal antibody against human Coagulation factor III/Tissue Factor (CD142, TF). Validated for WB, IF/ICC, and ELISA. Detects endogenous levels in human, mouse, and rat samples. Suitable for hemostasis and thrombosis research.

[Learn More](#)

Parameter	Specification
Product Name	Coagulation factor III/Tissue Factor Rabbit pAb
Remarks/Alias	TF; TFA; CD142; Coagulation factor III/Tissue Factor
Species	Human
Gene ID	2152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-251 of human Coagulation factor III/Tissue Factor (NP_001984.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13726
Protein Weight	33kDa
Shipping	Ice bag