

## Anti-Coagulation Factor X (F10) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P00742

Item Number: CM0015514



### Introduction

Rabbit polyclonal antibody targeting human coagulation factor X (F10). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse. Suitable for hemostasis, thrombosis, and inflammation research.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | F10 Rabbit pAb                                                                                                                          |
| Remarks/Alias    | FX; FXA; F10                                                                                                                            |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 2159                                                                                                                                    |
| GeneID           | 2159                                                                                                                                    |
| Immunogen        | Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 235-488 of human F10 (NP_000495.1). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                                                       |
| Cross Reactivity | Human, Mouse                                                                                                                            |
| SWISS            | P00742                                                                                                                                  |
| Protein Weight   | 55kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |