


Anti-Coch Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - O43405

Item Number: CM0030670



Introduction

Rabbit polyclonal antibody targeting human cochlin (COCH/COCH5B2), validated for IHC-P, IF/ICC, and ELISA. Expressed in the inner ear, cochlin regulates cell shape and motility. Ideal for hearing research.

[Learn More](#)

Parameter	Value
Product Name	COCH Rabbit pAb
Remarks/Alias	DFNA9; COCH5B2; DFNB110; COCH-5B2; COCH
Species	Human
Gene ID (human)	1690
Gene ID	1690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-260 of human COCH (NP_001128530.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O43405
Protein Weight	59kDa
Shipping	Ice bag