

Anti-Cofilin Rabbit Monoclonal Antibody For Wb, Elisa - P23528

Item Number: CM0008660



Introduction

CamelBio's rabbit anti-Cofilin monoclonal antibody targets human Cofilin-1 (P23528) with cross-reactivity to mouse and rat. Validated for WB and ELISA, it aids in studying actin dynamics, cell migration, and embryonic development.

[Learn More](#)

Parameter	Value
Product Name	Cofilin Rabbit mAb
Remarks/Alias	CFL; cofilin; HEL-S-15; Cofilin
Species	Human
Gene ID	1072 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cofilin (P23528)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23528
Protein Weight	19kDa
Shipping	Ice bag