

Anti-Cog1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wtw3

Item Number: CM0022327

Introduction

Rabbit polyclonal antibody against human COG1, a conserved oligomeric Golgi complex subunit. Suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, and rat. Ideal for Golgi function and CDG2G research.

[Learn More](#)

Parameter	Value
Product Name	COG1 Rabbit pAb
Remarks/Alias	LDLB; CDG2G; COG1
Species	Human
GeneID (Human)	9382
Gene ID	9382
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 835-935 of human COG1 (NP_061184.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTW3
Protein Weight	109kDa
Shipping	Ice bag