


Anti-Cog4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H9E3

Item Number: CM0014176



Introduction

Rabbit polyclonal antibody targeting human COG4, validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for studying Golgi function, SNARE-pin assembly, and retrograde transport.

[Learn More](#)

Parameter	Value
Product Name	COG4 Rabbit pAb
Remarks/Alias	COD1; CDG2J; SWILS; COG4
Species	Human
Gene ID (Human)	25839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-410 of human COG4 (NP_056201.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H9E3
Protein Weight	89kDa
Shipping	Ice bag