

Anti-Cog6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2V7

Item Number: CM0015469

Introduction

High-quality rabbit polyclonal antibody against human COG6 (SWISS Q9Y2V7), validated for Western blot and ELISA. Ideal for Golgi apparatus research and congenital glycosylation disorder studies.

[Learn More](#)

Parameter	Value
Product Name	COG6 Rabbit pAb
Remarks/Alias	COD2; SHNS; CDG2L; COG6
Species	Human
GenelD (human)	57511
GenelD	57511
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 508-657 of human COG6 (NP_065802.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2V7
Protein Weight	73kDa
Shipping	Ice bag