

Anti-Cog7 Rabbit Polyclonal Antibody For If/Icc, Elisa - P83436

Item Number: CM0028603

Introduction

Rabbit polyclonal antibody targeting human COG7, validated for IF/ICC and ELISA. Suitable for Golgi function and protein transport studies. Recognizes human COG7 (P83436).

[Learn More](#)

Parameter	Value
Product Name	COG7 Rabbit pAb
Remarks/Alias	CDG2E; COG7
Species	Human
GeneID (human)	91949
GeneID	91949
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 511-770 of human COG7 (NP_705831.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P83436
Protein Weight	86kDa
Shipping	Ice bag