

Anti-Coilin Polyclonal Antibody For Wb, Ip, Elisa - P38432

Item Number: CM0029420

Introduction

Rabbit polyclonal antibody against human Coilin (p80-coilin) for Western blot, immunoprecipitation, and ELISA. Cross-reacts with mouse and rat. Suitable for Cajal body dynamics and snRNP assembly studies.

[Learn More](#)

Parameter	Value
Product Name	Coilin Rabbit pAb
Remarks/Alias	CLN80; p80-coilin; Coilin
Species	Human
GeneID (Human)	8161
GeneID	8161
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 297-576 of human Coilin (NP_004636.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38432
Protein Weight	63kDa
Shipping	Ice bag