

Anti-Coilin Rabbit Monoclonal Antibody For Wb, Elisa - P38432

Item Number: CM0011882

Introduction

Recombinant rabbit monoclonal antibody against human Coilin (p80-coilin), validated for WB and ELISA. Detects Coilin in Cajal bodies; cross-reacts with mouse. Ideal for snRNP biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Coilin Rabbit mAb
Remarks/Alias	CLN80; p80-coilin; Coilin
Species	Human
GeneID (Human)	8161
GeneID	8161
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Coilin (P38432).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P38432
Protein Weight	63kDa
Shipping	Ice bag