


## Anti-Col10A1 (Collagen Alpha-1(X) Chain) Rabbit Polyclonal Antibody For Wb, Elisa - Q03692

Item Number: CM0022846



### Introduction

Rabbit polyclonal antibody against human COL10A1 (collagen alpha-1(X) chain), suitable for WB and ELISA. Useful for research on cartilage mineralization, endochondral ossification, and bone growth disorders.

[Learn More](#)

| Parameter        | Value                              |
|------------------|------------------------------------|
| Product Name     | COL10A1 Rabbit pAb                 |
| Remarks/Alias    | COL10A1; collagen alpha-1(X) chain |
| Species          | Human                              |
| GeneID (Human)   | 1300                               |
| GeneID           | 1300                               |
| Immunogen        | Recombinant protein                |
| Source           | Rabbit                             |
| Category         | Polyclonal Antibodies              |
| Application      | WB, ELISA                          |
| Cross Reactivity | Human                              |
| SWISS            | Q03692                             |
| Protein Weight   | 66kDa                              |
| Shipping         | Ice bag                            |