

Anti-Col11A1 Monoclonal Antibody For Wb, If/Icc, Elisa - P12107

Item Number: CM0005622

Introduction

Rabbit monoclonal antibody against human COL11A1, suitable for WB, IF/ICC, and ELISA. Ideal for extracellular matrix and collagen fibrillogenesis research. Sourced from CamelBio's one-stop IVD raw material supply.

[Learn More](#)

Parameter	Value
Product Name	COL11A1 Rabbit mAb
Remarks/Alias	STL2, COLL6, CO11A1, DFNA37, COL11A1
Species	Human
GeneID (Human)	1301
GeneID	1301
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 301-400 of human COL11A1 (NP_001845.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P12107
Protein Weight	167kDa/176kDa/181kDa/182kDa
Shipping	Ice bag