

Anti-Col11A1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P12107

Item Number: CM0013792

Introduction

Rabbit polyclonal antibody against human COL11A1, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for studying cartilage extracellular matrix and collagen fibrillogenesis.

[Learn More](#)

Parameter	Value
Product Name	COL11A1 Rabbit pAb
Remarks/Alias	STL2, COLL6, CO11A1, DFNA37, COL11A1
Species	Human
GeneID (Human)	1301
GeneID	1301
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-430 of human COL11A1 (NP_542196.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12107
Protein Weight	181kDa
Shipping	Ice bag