

Anti-Col11A2 Polyclonal Antibody For Wb, Elisa - P13942

Item Number: CM0013643

Introduction

COL11A2 Rabbit polyclonal antibody validated for WB and ELISA. Targets human collagen alpha-2(XI) chain (172 kDa), involved in fibrillogenesis. Useful for ECM research, hearing loss studies, and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	COL11A2 Rabbit pAb
Remarks/Alias	HKES, PARP, STL3, FBCG2, DFNA13, DFNB53, OSMEDA, OSMEDB, COL11A2
Species	Human
Gene ID (Human)	1302
Gene ID	1302
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-380 of human COL11A2 (NP_542411.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13942
Protein Weight	172kDa
Shipping	Ice bag