

Anti-Col1A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P02452

Item Number: CM0025577



Introduction

Rabbit polyclonal antibody against human COL1A1 (Collagen alpha-1(I) chain). Validated for WB, IF/ICC, and ELISA. Detects human and mouse COL1A1. Recombinant immunogen (aa 1359-1464). For extracellular matrix, collagen, and bone research.

[Learn More](#)

Parameter	Value
Product Name	COL1A1 Rabbit pAb
Remarks/Alias	OI1; OI2; OI3; OI4; EDSC; CAFYD; EDSARTH1; Collagen I/COL1A1
Species	Human
GeneID (Human)	1277
GeneID	1277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1359-1464 of human COL1A1 (NP_000079.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02452
Protein Weight	139kDa
Shipping	Ice bag