


Anti-Col2A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P02458

Item Number: CM0009194



Introduction

Validated COL2A1 Rabbit Monoclonal Antibody for WB, IHC-P, ELISA. Reacts with human COL2A1; essential for cartilage research and skeletal development studies.

[Learn More](#)

Parameter	Value
Product Name	COL2A1 Rabbit mAb
Remarks/Alias	AOM; ACG2; ANFH; LCPD; SEDC; STL1; VPED; ANFH1; EDMMD; OSCDP; PLSDT; SEDSTN; SMDALG; COL11A3; SEMDSTWK
Species	Human
GeneID (Human)	1280
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1242-1487 of human COL2A1 (NP_001835.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P02458
Protein Weight	30 kDa/134 kDa/142 kDa
Shipping	Ice bag