

Anti-Col2A1 Rabbit Polyclonal Antibody For Wb, Elisa - P28481

Item Number: CM0026642



Introduction

Rabbit polyclonal antibody against COL2A1, specific for cartilaginous tissues. Validated for WB and ELISA, reacts with mouse and rat. Ideal for skeletal development and cartilage research.

[Learn More](#)

Parameter	Details
Product Name	COL2A1 Rabbit pAb
Remarks/Alias	Dmm; Lpk; Col2; Del1; Col2a; Col2a-1; M100413; Rgsc413; Rgsc856
Species	Mouse
GeneID (Human)	1280
GeneID	12824
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P28481
Protein Weight	132-142 kDa
Shipping	Ice bag