


Anti-Col4A2 Polyclonal Antibody For Wb, Elisa - P08572

Item Number: CM0028943



Introduction

Rabbit polyclonal antibody targeting human collagen alpha-2(IV) chain (COL4A2). Validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Useful for basement membrane, angiogenesis, and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	COL4A2 Rabbit pAb
Remarks/Alias	ICH; BSVD2; POREN2; COL4A2
Species	Human
Gene ID (Human)	1284
Gene ID	1284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1486-1712 of human COL4A2 (NP_001837.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08572
Protein Weight	168kDa
Shipping	Ice bag