


Anti-Col4A3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q01955

Item Number: CM0021639



Introduction

COL4A3 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets collagen alpha-3(IV) chain, a key basement membrane protein involved in glomerular filtration and anti-angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	COL4A3 Rabbit pAb
Remarks/Alias	ATS2; ATS3; BFH2; COL4A3
Species	Human
GeneID (Human)	1285
GeneID	1285
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1555-1669 of human COL4A3 (NP_000082.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01955
Protein Weight	162kDa
Shipping	Ice bag