

Anti-Col4A3Bp Polyclonal Antibody (Cert, Gbp, Stard11) - Q9Y5P4

Item Number: CM0014668



Introduction

Rabbit polyclonal antibody against human COL4A3BP (CERT, GPBP, STARD11). Suitable for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Targets ceramide transfer protein mediating non-vesicular ER-to-Golgi lipid trafficking.

[Learn More](#)

Parameter	Value
Product Name	COL4A3BP Rabbit pAb
Remarks/Alias	CERT; GPBP; CERTL; MRD34; STARD11; COL4A3BP
Species	Human
Gene ID (Human)	10087
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 369-598 of human COL4A3BP (NP_112729.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5P4
Protein Weight	71kDa
Shipping	Ice bag