


Anti-Col4A4 Polyclonal Antibody For Wb And Elisa - P53420

Item Number: CM0022788



Introduction

Rabbit polyclonal antibody against human collagen alpha-4(IV) chain (COL4A4), suitable for WB and ELISA. Ideal for studying basement membrane biology, Alport syndrome, and glomerular function.

[Learn More](#)

Parameter	Value
Product Name	COL4A4 Rabbit pAb
Remarks/Alias	BFH; ATS2; BFH1; CA44; COL4A4
Species	Human
GeneID (Human)	1286
GeneID	1286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1573-1690 of human COL4A4. (NP_000083.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P53420
Protein Weight	164kDa
Shipping	Ice bag