


Anti-Col4A5 Polyclonal Antibody For Wb, If/Icc, Elisa - P29400

Item Number: CM0028969



Introduction

Rabbit polyclonal antibody against human type IV collagen alpha-5 chain, validated for Western blot, IF/ICC, and ELISA. Ideal for investigating glomerular basement membrane biology and Alport syndrome.

[Learn More](#)

Parameter	Value
Product Name	COL4A5 Rabbit pAb
Remarks/Alias	ATS, ASLN, ATS1, CA54, COL4A5
Species	Human
Gene ID (Human)	1287
Gene ID	1287
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1462-1691 of human COL4A5 (NP_203699.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P29400
Protein Weight	161kDa
Shipping	Ice bag