

Anti-Col6A2 Polyclonal Antibody For Wb, Elisa - P12110

Item Number: CM0024702



Introduction

High-quality rabbit polyclonal antibody targeting human COL6A2. Validated for WB and ELISA, cross-reacts with mouse. Ideal for extracellular matrix and collagen VI research. SWISS: P12110.

[Learn More](#)

Parameter	Value
Product Name	COL6A2 Rabbit pAb
Remarks/Alias	UCMD1; BTHLM1; PP3610; COL6A2
Species	Human
GeneID (Human)	1292
GeneID	1292
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 720-1019 of human COL6A2 (NP_001840.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P12110
Protein Weight	109kDa
Shipping	Ice bag