

Anti-Col9A1 Polyclonal Antibody For Wb And Elisa - P20849

Item Number: CM0029357



Introduction

Rabbit polyclonal antibody against human COL9A1, validated for Western blot and ELISA with cross-reactivity to rat. Suitable for cartilage and vitreous extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	COL9A1 Rabbit pAb
Remarks/Alias	MED; EDM6; STL4; DJ149L1.1.2; COL9A1
Species	Human
Gene ID	1297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-270 of human COL9A1 (NP_001842.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P20849
Protein Weight	92kDa
Shipping	Ice bag