


Anti-Colec11 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bwp8

Item Number: CM0029650



Introduction

COLEC11 Rabbit Polyclonal Antibody suitable for WB, IHC-P, ELISA applications in human, mouse, rat. Targets collectin-11, a lectin involved in innate immunity, apoptosis, and embryogenesis. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	COLEC11 Rabbit pAb
Remarks/Alias	3MC2; CLK1; CL-11; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; COLEC11
Species	Human
Gene ID (Human)	78989
Gene ID	78989
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-271 of human COLEC11 (NP_076932.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BWP8
Protein Weight	29kDa
Shipping	Ice bag