

Anti-Colec12 Polyclonal Antibody For Wb, Elisa - Q5Ku26

Item Number: CM0013737



Introduction

COLEC12 Rabbit polyclonal antibody suitable for Western blot and ELISA. Reacts with human and mouse COLEC12 (CL-P1, SRCL). Useful for studying scavenger receptor functions in host defense and oxLDL uptake.

[Learn More](#)

Parameters	Details
Product Name	COLEC12 Rabbit pAb
Alias	CLP1, NSR2, SRCL, SCARA4, COLEC12
Species	Human
Gene ID (Human)	81035
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-270 of human COLEC12 (NP_569057.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5KU26
Protein Weight	82kDa
Shipping	Ice bag