

Anti-Colgalt2 Polyclonal Antibody For Wb, Elisa - Q8lyk4

Item Number: CM0014035

Introduction

Rabbit polyclonal antibody targeting procollagen galactosyltransferase 2 (COLGALT2). Validated for WB and ELISA with human and mouse cross-reactivity. Ideal for collagen glycosylation studies.

[Learn More](#)

Parameter	Value
Product Name	COLGALT2 Rabbit pAb
Remarks/Alias	C1orf17; GLT25D2; ColGalT 2; COLGALT2
Species	Human
Gene ID (Human)	23127
Gene ID	23127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-230 of human COLGALT2 (NP_055916.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IYK4
Protein Weight	73kDa
Shipping	Ice bag