

Anti-Collagen Alpha-1(Xvi) Chain (Col16A1) Rabbit Polyclonal Antibody For Wb And Elisa - Q07092

Item Number: CM0021968

Introduction

COL16A1 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat collagen alpha-1(XVI) chain, a secreted extracellular matrix protein involved in cell adhesion and integrin signaling.

[Learn More](#)

Parameter	Value
Product Name	COL16A1 Rabbit pAb
Remarks/Alias	447AA; FP1572; COL16A1
Species	Human
GeneID	1307
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1000-1100 of human COL16A1 (NP_001847.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07092
Protein Weight	158kDa
Shipping	Ice bag