


Anti-Collagen I/Col1A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P02452

Item Number: CM0008402



Introduction

High-quality rabbit monoclonal antibody against Collagen I/COL1A1, suitable for WB, IF-P, IHC-P, and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for extracellular matrix and bone biology research.

[Learn More](#)

Parameter	Value
Product Name	Collagen I/COL1A1 Rabbit mAb
Remarks/Alias	OI1; OI2; OI3; OI4; EDSC; CAFYD; EDSARTH1
Species	Human
Gene ID (Human)	1277
Gene ID	1277
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02452
Protein Weight	139kDa
Shipping	Ice bag