

Anti-Collagen I/Col1A1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P02452

Item Number: CM0018116



Introduction

Rabbit polyclonal antibody targeting human Collagen I/COL1A1, validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat tissues and cells. Suitable for extracellular matrix, bone and connective tissue studies.

[Learn More](#)

Parameter	Value
Product Name	Collagen I/COL1A1 Rabbit pAb
Remarks/Alias	OI1; OI2; OI3; OI4; EDSC; CAFYD; EDSARTH1; Collagen I/COL1A1
Species	Human
Gene ID (Human)	1277
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1151-1250 of human Collagen I/COL1A1 (NP_000079.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02452
Protein Weight	139 kDa
Shipping	Ice bag