

Anti-Collagen I/Col1A1 Rabbit Polymab® For Wb - P02452

Item Number: CM0006868



Introduction

Rabbit polymAb® against Collagen I/COL1A1 for WB, IHC-P, and IF applications. Detects human, mouse, and rat type I collagen with 139 kDa molecular weight. Ideal for extracellular matrix studies.

[Learn More](#)

Parameter	Value
Product Name	Collagen I/COL1A1 Rabbit PolymAb®
Remarks/Alias	OI1; OI2; OI3; OI4; EDSC; CAFYD; EDSARTH1
Species	Human
GeneID (human)	1277
GeneID	1277
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1150-1250 of human Collagen I/COL1A1 (NP_000079.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02452
Protein Weight	139 kDa
Shipping	Ice bag