

Anti-Collagen I/Col1A2 Monoclonal Antibody For Wb, Elisa - P08123

Item Number: CM0002420



Introduction

Collagen I/COL1A2 Rabbit mAb for WB, ELISA applications. Validated in human, mouse, rat. Target molecular weight 129 kDa. Ideal for extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	Collagen I/COL1A2 Rabbit mAb
Remarks/Alias	OI4; EDSCV; EDSARTH2; Collagen I/COL1A2
Species	Human
GeneID (Human)	1278
GeneID	1278
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500-600 of human Collagen I/COL1A2 (NP_000080.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08123
Protein Weight	129kDa
Shipping	Ice bag