

Anti-Collagen I/Col1A2 Monoclonal Antibody For Wb, If/Icc, Elisa - P08123

Item Number: CM0006551



Introduction

Rabbit monoclonal antibody targeting human Collagen I/COL1A2, validated for WB, IF/ICC, and ELISA. Suitable for studying extracellular matrix, tendon, ligament, and bone biology. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	Collagen I/COL1A2 Rabbit mAb
Remarks/Alias	O14; EDSVCV; EDSARTH2; Collagen I/COL1A2
Species	Human
GeneID (Human)	1278
GeneID	1278
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08123
Protein Weight	129kDa
Shipping	Ice bag