

Anti-Collagen I/Col1A2 Polyclonal Antibody For Wb, If/Icc, Elisa - P08123

Item Number: CM0021019



Introduction

Rabbit polyclonal antibody targeting human Collagen I/COL1A2, suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Valuable for extracellular matrix, bone, and connective tissue research.

[Learn More](#)

Parameter	Value
Product Name	Collagen I/COL1A2 Rabbit pAb
Remarks/Alias	OI4; EDSCV; EDSARTH2; Collagen I/COL1A2
Species	Human
GeneID	1278 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1161-1260 of human Collagen I/COL1A2 (NP_000080.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08123
Protein Weight	129kDa
Shipping	Ice bag