

Anti-Collagen I/Col1A2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P08123

Item Number: CM0027501



Introduction

High-quality rabbit polyclonal antibody against human Collagen I/COL1A2, validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying collagen biology, extracellular matrix, and bone mineralization.

[Learn More](#)

Parameter	Details
Product Name	Collagen I/COL1A2 Rabbit pAb
Remarks/Alias	OI4; EDSCV; EDSARTH2; Collagen I/COL1A2
Species	Human
GeneID (Human)	1278
GeneID	1278
Immunogen	Synthetic peptide corresponding to amino acids 501-600 of human Collagen I/COL1A2 (NP_000080.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08123
Protein Weight	129 kDa
Shipping	Ice bag