


Anti-Collagen II Polyclonal Antibody For Wb, Elisa - P02458

Item Number: CM0023105



Introduction

High-quality rabbit polyclonal antibody against human Collagen II (COL2A1). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for cartilage development and skeletal biology studies.

[Learn More](#)

Parameter	Value
Product Name	Collagen II Rabbit pAb
Remarks/Alias	AOM; ANFH; SEDC; STL1; COL11A3; Collagen II
Species	Human
Gene ID (Human)	1280
Gene ID	1280
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human Collagen II (NP_001835.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02458
Protein Weight	142kDa
Shipping	Ice bag