

Anti-Collagen Iii Alpha 1/Col3A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P02461

Item Number: CM0024992



Introduction

Rabbit polyclonal antibody targeting human Collagen III alpha 1 (COL3A1). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for extracellular matrix and connective tissue research.

[Learn More](#)

Parameter	Specification
Product Name	Collagen III alpha 1/COL3A1 Rabbit pAb
Remarks/Alias	EDS4A; EDSVASC; PMGEDSV; Collagen III alpha 1/COL3A1
Species	Human
GeneID (human)	1281
GeneID	1281
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1201-1301 of human Collagen III alpha 1/COL3A1 (NP_000081.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02461
Protein Weight	139kDa
Shipping	Ice bag