

Anti-Collagen Iii Alpha 1/Col3A1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P02461

Item Number: CM0003307



Introduction

High-quality rabbit monoclonal antibody against Collagen III alpha 1/COL3A1 (P02461). Validated for western blot, immunoprecipitation, and ELISA. Reacts with human, mouse, and rat. A reliable tool for extracellular matrix and connective tissue research.

[Learn More](#)

Parameter	Value
Product Name	Collagen III alpha 1/COL3A1 Rabbit mAb
Remarks/Alias	EDS4A; EDSVASC; PMGEDSV; Collagen III alpha 1/COL3A1
Species	Human
Gene ID (Human)	1281
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-400 of human Collagen III alpha 1/COL3A1 (P02461)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02461
Protein Weight	139kDa
Shipping	Ice bag