


Anti-Collagen Iv Monoclonal Antibody For Wb, Db, Elisa - P53420 / Q14031 / P08572 / P02462

Item Number: CM0006635



Introduction

Collagen IV Rabbit Monoclonal Antibody validated for Western blot, dot blot, ELISA. Cross-reacts with Human, Mouse, Rat. Useful for basement membrane and glomerular disease research.

[Learn More](#)

Parameter	Value
Product Name	Collagen IV Rabbit mAb
Remarks/Alias	BFH; ATS2; BFH1; CA44; Collagen IV
Species	Human
Gene ID (Human)	1286/1288/1284/1282
Gene ID	1286/1288/1284/1282
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1573-1690 of human COL4A4 (NP_000083.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53420/Q14031/P08572/P02462
Protein Weight	164 kDa
Shipping	Ice bag