


Anti-Collagen Iv Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q01955 / P29400 / P08572 / P02462

Item Number: CM0004951



Introduction

Rabbit monoclonal antibody targeting Collagen IV alpha-3 chain (Goodpasture antigen). Validated for WB, IF/ICC, ELISA in human and rat samples. Suitable for basement membrane and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Collagen IV Rabbit mAb
Remarks/Alias	ATS2; ATS3; BFH2; COL4A3/COL4A5/COL4A2/COL4A1
Species	Human
Gene ID (Human)	1285/1287/1284/1282
Gene ID	1285/1287/1284/1282
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1555-1669 of human COL4A3 (NP_000082.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q01955/P29400/P08572/P02462
Protein Weight	162kDa
Shipping	Ice bag