


Anti-Collagen VI/Col6A1 Polyclonal Antibody For Wb, Elisa - P12109

Item Number: CM0028234



Introduction

Rabbit polyclonal antibody against human Collagen VI/COL6A1, validated for WB and ELISA applications. Cross-reacts with mouse. Ideal for studies of extracellular matrix, cell adhesion, and muscular dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	Collagen VI/COL6A1 Rabbit pAb
Remarks/Alias	OPLL; BTHLM1; UCHMD1; Collagen VI/COL6A1
Species	Human
GeneID (Human)	1291
GeneID	1291
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-250 of human Collagen VI/Collagen VI/COL6A1 (NP_001839.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P12109
Protein Weight	109kDa
Shipping	Ice bag